



# Meta-Analysis of 33 Randomized Controlled Trials Demonstrates Value of CAde in Colorectal Polyp Detection

## Research Goal

Evaluate the potential advantages and disadvantages of use of computer-aided detection (CAde) during colorectal cancer (CRC) screening colonoscopy.

Systematic review and meta-analysis

33 randomized controlled trials (RCTs)

## Key Results

CAde significantly improved detection of :

- Advanced adenomas
- Diminutive and small adenomas
- Proximal and distal adenomas

## Clinical Relevance

- Colonoscopy is the gold standard of CRC screening, allowing for the identification and removal of pre-cancerous lesions before they progress.<sup>2</sup>
- While >10mm adenomas have the greatest potential for malignancy, diminutive (0-5mm) and small lesions (6-10mm) may also develop into CRC.
  - 0.6-5.6% of <10mm polyps show high-grade dysplasia and malignancy.
  - 3+ non-advanced diminutive adenomas have been associated with elevated metachronous cancer risk.
- Interval CRCs are more common in the proximal colon.

↑ 1% increase in ADR

↓ 3% decrease in risk of interval CRC

“Our study confirmed the role of AI-aided colonoscopy in safely improving the detection of colorectal neoplasia...without significant delays in procedure time.”

– Lou et. al 2023

Artificial intelligence for colorectal neoplasia detection during colonoscopy: A systematic review and meta-analysis of randomized clinical trials<sup>1</sup>

Shenghan Lou, Fenqui Du, Wenjie Song, Yixiu Xia, Xinyu Yue, Da Yang, Binbin Cui, Yanlong Liu, and Peng Han Originally published in eClinical Medicine Volume 66 (December 2023)

ADR = adenoma detection rate



Scan to read the study

# Methodology

Meta-analysis of RCTs comparing CADe to standard of care

## Primary outcome measures:

adenoma miss rate (AMR), adenoma detection rate (ADR), and adenomas per colonoscopy (APC)

## Secondary outcome measures:

polyp miss rate (PMR), polyp detection rate (PDR), polyps detected per colonoscopy (PPC), procedure time, adverse events, and false alarms



**33 unique  
randomized  
controlled trials**

Including 8 tandem studies



**27,404 patients**

Excluded trials of patients with inflammatory bowel disease or hereditary polyposis syndromes

## Results

AI significantly **decreased polyp and adenoma miss rates**

**↓ 51%** relative decrease  
in **AMR**

**↓ 53%** relative decrease  
in **PMR**

AI significantly **increased polyp and adenoma detection rates**

(1 additional adenoma for every ~5 patients examined)

**↑ 24%** relative increase  
in **ADR & PDR**

**↑ 39%** relative increase  
in **APC & PPC**

“Concerns regarding increased procedure time and complication risk with AI-aided colonoscopy were not substantially supported by our findings... adding an AI-aided system will not significantly affect the colonoscopy workflow”

– Lou et. al 2023